

Delta-Lite® Plus

Fiberglass Cast Tape



**The Trusted
Rigid Solution**

Comfort and Ease of Application

Should Go Without Saying

Delta-Lite® Plus is designed to ensure excellent rigid support. Thanks to its light but strong fiberglass substrate and resin formula it has become the trusted solution for all fractures requiring strong and lasting support.


Strong yet light fiberglass yarns

- Rigid and durable casts for uncompromising protection
- Substrate is still flexible enough to avoid muscle constrictions
- No premature cast changes due to cast break down

Excellent lamination for lasting support

- Secure lamination between the individual layers for long-lasting stability and support
- Weight-bearing in approximately 20 minutes
- Creates long-lasting, lightweight cast
- Noticeable stability conveys trust

Secure Lamination



Strong Support

Delta-Lite® Plus

Fiberglass Cast Tape

While strong support and durability are essential features of a cast, wearing comfort and ease of application are even more apparent to patient and applicator. Therefore, Delta-Lite® Plus comes with a number of features to enhance these properties.

Resin formula with slip agent surfactant

- Slip agent surfactant creates a plaster of paris-like sensation during working time for precise anatomic molding
- Reliable, comfortable working time window
- Good pliability throughout working time
- Freedom to use glove of your choice

Smooth finish improves wearing comfort

- Avoids snag on clothes and bed linens
- Excellent end lay-down for even, smooth finish
- Improves patient comfort and compliance



Delta-Lite® Plus

Fiberglass Cast Tape

Proven Overall Quality

In a recent study among 387 clinicians in the USA, Delta-Lite® Plus was compared to the cast tape currently used by each clinician. Nine key product features, such as **lamination**, **rigidity**, **conformability** and **smooth finish** were chosen for comparison. In all nine categories Delta-Lite® Plus was rated higher than the currently used fiberglass tapes. **Over 90% of all clinicians in the study rated Delta-Lite® Plus' overall quality superior to their currently used product.**

Code	Dimensions	Box Contents
□ 7345800	2.5 cm x 1.8 m	10 rolls
□ 7345801	5 cm x 3.6 m	10 rolls
□ 7345802	7.5 cm x 3.6 m	10 rolls
□ 7345803	10 cm x 3.6 m	10 rolls
□ 7345804	12.5 cm x 3.6 m	10 rolls
■ 7345820	5 cm x 3.6 m	10 rolls
■ 7345821	7.5 cm x 3.6 m	10 rolls
■ 7345822	10 cm x 3.6 m	10 rolls
■ 7345823	12.5 cm x 3.6 m	10 rolls
■ 7345825	5 cm x 3.6 m	10 rolls
■ 7345826	7.5 cm x 3.6 m	10 rolls
■ 7345827	10 cm x 3.6 m	10 rolls
■ 7345828	12.5 cm x 3.6 m	10 rolls
■ 7345830	5 cm x 3.6 m	10 rolls
■ 7345831	5 cm x 3.6 m	10 rolls
■ 7345832	10 cm x 3.6 m	10 rolls
■ 7345835	5 cm x 3.6 m	10 rolls
■ 7345836	7.5 cm x 3.6 m	10 rolls
■ 7345837	10 cm x 3.6 m	10 rolls
■ 7345838	12.5 cm x 3.6 m	10 rolls

Code	Dimensions	Box Contents
■ 7345845	5 cm x 3.6 m	10 rolls
■ 7345846	7.5 cm x 3.6 m	10 rolls
■ 7345847	10 cm x 3.6 m	10 rolls
□ 7345850*	5 cm x 3.6 m	10 rolls
□ 7345851*	7.5 cm x 3.6 m	10 rolls
□ 7345852*	10 cm x 3.6 m	10 rolls
■ 7345855	5 cm x 3.6 m	10 rolls
■ 7345856	7.5 cm x 3.6 m	10 rolls
■ 7345857	10 cm x 3.6 m	10 rolls
■ 7345860	5 cm x 3.6 m	10 rolls
■ 7345861	7.5 cm x 3.6 m	10 rolls
■ 7345862	10 cm x 3.6 m	10 rolls
■ 7345865	5 cm x 3.6 m	10 rolls
■ 7345866	7.5 cm x 3.6 m	10 rolls
■ 7345867	10 cm x 3.6 m	10 rolls
■ 7345875	5 cm x 3.6 m	10 rolls
■ 7345876	7.5 cm x 3.6 m	10 rolls
■ 7345877	10 cm x 3.6 m	10 rolls

Colors and availability subject to change. Please refer to the product label or package insert for full instructions on the safe use of these products.

*Mixed pack contains 2 rolls of each: deep blue, red, light blue, purple, pink.

In combination with Delta-Dry® water resistant cast padding, Delta-Lite® Plus can be upgraded to a fully water-resistant cast, allowing patients to bathe and shower.

