

# INRA Freeze®

## The first equine semen freezing extender

**INRA Freeze®:**  
**increase fertility results!**

Ease your daily work - ready to use

Transport at room temperature

Avoid sanitary risks - sterile

**Ready to use**

**NEW**

**INRA Freeze®:**  
**improve your freezing techniques!**

- Better semen protection due to the association of INRA 96 and egg yolk plasma
- Proven to increase fertility in vivo: results published by INRA
- Simple to use:  
INRA 96® as centrifugation extender  
+ INRA Freeze® as freezing extender



WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES




[www.imv-technologies.com](http://www.imv-technologies.com)

# INRA FREEZE®:

A significant fertility increase proven in vivo

Research completed by Institut National de la Recherche Agronomique (INRA) have shown a real **in vivo fertility performances increase**

Post-thaw, insemination with INRA Freeze® significantly increases performance and fertility compared to the traditional extender



PILLET E.\* et al., 2008.

Freezing stallion semen in INRA96® based extender improves fertility rates in comparison to INRA 82, Dairy Science and Technology, vol.88, n 2, p257-265.

PILLET E.\* et al., 2009.

A new freezing extender for stallion semen to get high fertility rate: INRA Freeze.

Reward: Young Researcher, 60<sup>th</sup> Annual Meeting of the European Federation of Animal Science, Barcelona, August 2009.

**Protocol :** 154 cycles inseminated with 8 straws ( x 0.5 ml ) . AI dose = 400 million total sperm.

Motility > 35% post thaw.

**One AI at 30 hours after induced ovulation**

## INRAFreeze®, the unique equine semen freezing extender

- Give better fertility results
- Easy to use and ready to use
- Easy to transport
- Easy to store

### Equine Freezing Pack

Ref 022110

Includes 2 INRA Freeze (400 ml), 2 INRA 96 (240 ml) and 1 MAXIFreeze (20 ml)

WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES



www.imv-technologies.com